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
• Protocol: microplate

1. Pipette 25 µl of each standard or unknown sample replicate into a microplate well (working range = 20-2,000µg/ml).

Note: If sample size is limited,10µl of each unknown sample and standard can be used sample to Working Solution; WS ratio=1:20).However, the working range of the assay in this case will be limited to 125-2,000µg/ml.

2. Add 200 µl of the Working Solution (WS) to each well and mix plate thoroughly on a plate shaker for 30seconds.

Note: To preparation of Working Solution(WS), see the "PREPARATION of STANDARDS and WORKING SOLUTION"

3. Cover plate and incubate at 37°C for 30minutes.
4. Cool plate to room temperature(RT).
5. Measure the absorbance at or near 562 nm on a plate reader.


Note: Wavelengths from 540-590nm have been used successfully with this method. If higher 562nm measurements are desired, increase the incubation time to 2-hours.

6. Subtract the average 562 nm absorbance measurement of the Blank standard replicates from the 562nm measurement so fall other individual standard and unknown sample replicates

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
EXPERIMENTAL DATA

IMPROVED SENSITIVITY & LINEAR WORKING RANGE

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